

new mom™

ELASTIC ABDOMINAL CORSET 4 PANEL



- Stretchable material allows easy customization as the waist size reduces

Intended Use of the Product

To provide adequate compression and support to the abdominal region post surgery/delivery

Directions for use

Wrap the product around the abdomen

Stretch out both the ends of the product and fasten the hook & loop

Warning

If discomfort, pain or irritation results, persists or increases, discontinue use and consult physician

Caution⚠

To be used under the advice and supervision of a registered medical practitioner

Dispose after use as per hospital policies and/ or country's regulatory norms

Product care



LATEX FREE

Size Available

Measure around the hip

One size fits most
70-130cm



For more
Products
Scan
QR Code



ISO 13485:2016
Certified company



Dynamic Techno Medicals Pvt. Ltd.
Kodikuthimala, Asokapuram, Aluva-683 101, Ernakulam, Kerala, India
Toll Free Number:- 1800 102 7902, www.dynamictechnomedicals.com
MADE IN INDIA